

# **Status Update On The New HMS Wire Chamber**

## **Joint Hall A & C Summer Collaboration Meeting**

### **June 22-23, 2017**

**Bishnu Pandey**  
**Hampton University, VA, USA**



# Outline

- Introduction
- Basic Design of HMS Chamber
- Test of HMS Chamber
- Cosmic Rays Signal Test Set Up
- Plateau Curve
- Hardware Status at ESB
- Test Summary at ESB
- Summary and Future Plans
- Some Snapshots
- Acknowledgements

# Introduction

The new HMS Chambers are being built to replace the old existing chambers using the same design of the SHMS chambers. The chambers will be used to track the scattered particles after beam collides with a fixed target to determine the momentum and scattering angle.

# Chamber Construction at Hampton University

## First Chamber

Started: May 10 2016

Completed: August 10 2016

Shipped to Jlab: 08/12/2016

## Second Chamber

Started: August 12 2016


completed: November 06 2016

Shipped to Jlab: 11/11/2016

# Basic Design

It consists of:



- 2 cathode windows (KW)
- 8 cathode planes (K)
- 6 wire planes ( X,X',U,U',V,V')
- 1 middle plane with card carriers and readout electronics
- 2 Aluminum frames



# Construction of Planes

## Cathode plane

foil is stretched along a table connected with a vacuum pump and a k plane after applying mechanical epoxy is placed over the foil for 24 hours.







# Wire plane

Wires are strung along the wire plane over the printed circuit board. Tension is fixed by hanging an iron bar on the wire. The position of the wire is fixed by solder and epoxy.

U-Plane:



X-Plane:




Picture shows SHMS wire plane

## Types of wire used

- 20  $\mu\text{m}$  Gold plated Tungsten sense wire
- 100  $\mu\text{m}$  Copper-Beryllium field wire



# Some Other Elements of Chamber





# Operational Testing at ESB

## Gas mixture used

$Ar : CO_2$  75:25 by volume

$Ar$  Ionizing gas


$CO_2$  quenching gas

The mixture is non flammable

# Conditioning

High voltage conditioning started before mid November and completed around first week of January 2017.


- chamber is connected to gas
- raised voltage up to 1860V
- At 1800V current per plane is less than 100nA





# Cosmic Ray Signal Test Set up

Started signal test beginning in January 2017.






# Plateau Curve


Plateau measurement started mid January 2017.

Plateau measurement with  $V_{th} = 3V$  (Chamber I)





# Plateau measurement with $V_{th} = 3V$ (Chamber II)

Gas: Ar/Co<sub>2</sub>




# Hardware Status at ESB

- All the ribbon cables connected from chamber (Nanometric cards) to the TDCs.
- Low voltage, threshold voltage and the high voltage are connected.
- Gas mixture,  $ArCO_2$  is connected to the chamber.






# Test of Chamber




# Test Summary at ESB

For threshold voltage 5.5 volt

## HDC 1X1 Wiremap




# HDC 1X2 Wiremap







# HDC 1U1 Wiremap




# HDC 1U2 Wiremap



# HDC 1V1 Wiremap



# HDC 1V2 Wiremap



# Summary


- Two HMS wire chambers were built at Hampton University from May to November 2016.
- After the primary test at Hampton, the chambers were shipped to Jlab (ESB building) for testing.
- High voltage conditioning completed mid-January 2017.
- Plateau curve test is done and plateau voltage found around 1750 volt.
- Started taking data from third week of June 2017.



# Future Plans

- The signal from each wire and the efficiency of each chamber will be tested.
- Chambers to be mounted on the frame in ESB and fiducialized.
- Old Drift chambers will be removed from HMS detector stack late July 2017.
- New Chambers will be installed in Hall C around mid August 2017.
- Cosmic testing in Hall C using  $ArCO_2$





calorimeter

Hodoscope


Cherenkov  
detector

Vacuum vessel

Old existing  
chamber




# Particle Detectors inside the HMS







# Some Snap shots









Detector  
frame

Ribbon cable  
with label

Ch2X'2(12-27)  
card#4

# Acknowledgements

Dr. Liguang Tang

Dr. Brad Sawatzky

Dr. Steve Wood

Mahlon Long (Chuck)

Thank You